

I-View Gold

CMOS INTRAORAL SENSOR

Trident proudly introduces the most advanced technology to obtain intraoral images in less time with less radiation exposure.

I-ViewGold is perfect to get better results in no time. The sensor covers a wide range of exposure setting, capturing a consistent and diagnosable image under any conditions.


Combine I-View Gold with your traditional X-Ray generator and obtain incredibly detailed images.

Two sizes

Thinner body and even surface

Rounded corners for a comfortable and easy positioning

Water-resistant



Size 1
20 x 30 mm
active area

Size 2
26 x 36 mm
active area

Practical cover to protect the active area from scratches and damages.



trident

I-View Gold



Dimensions

Size 1 25.4 x 36.8mm
Size 2.0 31.3 x 42.9mm



Specifications

Detector CMOS
Pixel size 18 μ m
Theoretical Resolution 28 lp/mm
Dynamic Range 12 bit



Deep View 4[®] Software

New generation of Deep-View software, designed to efficiently create, organize store and share, your patients' data files and images

- Multiple database management
- **Easy images acquisition device configuration**
- Advanced image filters
- **Wide choice and customization of templates for Full Mouth Series exams**
- Simulation of implant measurement and medical reports with teething preview
- **Voice guidance option for image acquisition**
- Images import/export in DICOM and other graphic formats
- **Printing images with layout customization**

Direct USB connection
2.7 m cable length
Tear-resistant and flexible sensor body's connection

trident

www.trident-dental.com



DENTALTI
| simplifying dental technology

1.800.672.5733
www.dentalti.com